

CPB Initiation

1

Closed loop communication

"Going on bypass"

2

Start timer

3

Communication

"Arterial on"

- Open arterial clamp
- Flow less than 500 mL/min

4

Communication

"Venous Open"

- Open venous clamp slowly
- Observe reservoir volume ↑ increasing, e.g. draining patient

5

Increase flow / RPMs + increase venous drainage from patient

See Flow Ramp — Worked Example

6

Confirm

- e2.4 cl
- Good drainage
- Zero CVP
- No ejection
- Good line pressure
- Colour difference in lines

7

Communication

"Full flow"

FLOW RAMP — WORKED EXAMPLE

Flow	Venous open	CVP	Arterial waveform	Line pressure
1/4	25%	8		80
1/2	50%	6		110
3/4	75%	4		150
Full	100%	0		180

Flat line = zero ejections

FULL-FLOW CONFIRMATION

- ☐ e2.4 cl
- ☐ Good drainage
- ☐ Zero CVP
- ☐ No ejections
- ☐ Good line pressure
- ☐ Colour difference in lines

THESE MUST BE CONFIRMED BEFORE COMMUNICATIONS